## Intended use

Hygicult E/β-GUR slides are intended for presumptive detection of bacteria belonging to the family *Enterobacteriaceae* and for identification of species producing β-glucuronidase enzyme. The slide is covered on one side with modified VRB Agar (Violet Red Bile Agar with addition of glucose) which allows bacteria belonging to *Enterobacteriaceae* to grow as red colonies. The glucose also allows some other species to grow as red colonies. The other side of the slide is covered with colourless β-Gur Agar. Species producing β-glucuronidase are able to grow on this side as brown colonies. The test can be performed on-site for monitoring different types of materials, both solid and liquid. As required, the slides can be used as convenient transport media for samples. **Note:** The limit values for microbial count in normal drinking water are too low to be detected by the Hygicult method.

## Contents of the kit

| Hygicult E/β-GUR     | Cat. No. 68267 |
|----------------------|----------------|
| Test slides          | 10 pcs         |
| Labels               | 10 pcs         |
| Instructions for use | 1 pc           |

## Typical formulation

| Modified VRB Agar | β-GUR Agar                     |
|-------------------|--------------------------------|
| Peptone           | Sodium lauryl sulphate         |
| Yeast extract     | 8-hydroxyquinoline glucuronide |
| Sodium chloride   | Yeast extract                  |
| Lactose           | Vitamins                       |
| Glucose D         | Cysteine hydrochloride         |
| Bile salts        | Magnesium sulphate             |
| Neutral red       | Bile salts                     |
| Crystal violet    | Manganese (II) chloride        |
| Agar agar         | Iron (III) citrate             |
| Water             | Agar agar                      |
|                   | Water                          |

## Warnings and precautions

Do not use product beyond the expiry date marked on the kit. Do not use the kit if you notice

- discoloration or dehydration of the growth medium
- detachment of the growth media from the plastic slide
- evidence of bacterial or fungal growth

Do not touch the growth because any colony growing on the slide may be pathogenic.

## Storage

Store the kit at room temperature (18...25°C / 64...77°F) protected from draught, temperature fluctuations and light sources. Avoid storage near heat-generating appliances. Do not allow to freeze. The expiry date (year-month-date) is marked on the box and on the cap of each slide.

## Sampling

To avoid contamination, the growth medium should not come into contact with any other material than the one to be tested. On the other hand, it is important that the growth medium makes full contact with the material to be tested. After sampling screw the slide tightly back into the tube.

## Contact inoculation (Fig. 1a, 1b)

Solid surfaces can be tested by pressing each side of the slide firmly against the surface for three or four seconds. The slide should be held still during pressing. The hinged design offers ease of use.

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## Dipping (Fig. 2)

Fluid samples are tested by dipping the slide in the liquid for three or four seconds. Blot the last drops on absorbent paper.

## Swabbing (Fig. 3)

Semisolid materials or objects that are difficult to reach can be tested by carefully rolling a sterile swab over an area delimited using e.g. a frame. If the object is dry, the swab should first be moistened with sterile water. The moistened swab can also be used for obtaining samples from powders (e.g. spices) or viscous fluids. After swabbing the sample area, roll the swab gently over the agar surfaces of the slide from left to right and from bottom to top.

## Incubation (Fig. 4)

Incubate the slide tightly enclosed in its tube at 35...37°C for 24–48 hours.

## Interpretation of results (Fig. 5)

Remove the slide from its tube after incubation and determine the microbial count (number of colony forming units, CFU) and examine the colour reactions by comparing with the model chart. Bacteria belonging to *Enterobacteriaceae* grow on the modified VRB Agar as red colonies. The glucose also allows some other gram-negative bacteria, e.g. *Pseudomonas* species, to grow as red colonies. β-glucuronidase-positive organisms grow on β-GUR agar as colonies in various shades of brown. As the colouration may be weak for some strains at high densities (10<sup>6-7</sup> CFU/ml), any shade of brown is indicative of β-glucuronidase-positive growth. β-glucuronidase activity is found in ca. 90 % of strains of *Escherichia coli*. Some species of *Salmonella*, *Edwardsiella*, *Shigella* and *Yersinia* are also β-glucuronidase producers. Gram-negative strains without β-glucuronidase activity grow as colourless colonies on this agar. The growth of gram-positive organisms is inhibited on both media. The following levels can be considered as a rough basis for evaluating the degree of contamination:

|                   | Contact inoculation |
|-------------------|---------------------|
| Clean             | 0 CFU/slide         |
| Contaminated      | 1–10 CFU/slide      |
| Very contaminated | > 10 CFU/slide      |

The presence of *Enterobacteriaceae* in cooked food always indicates mishandling of the product or inadequate hygiene.

## Limitations of the method

When used as a contact slide, Hygicult E/β-GUR equals the contact plate method in sensitivity, whereas the dip and swab procedures have a detection limit of 1000 CFU/ml. The allowed total microbial concentration of normal drinking water is too low to be reliably detected using Hygicult E/β-GUR. Results obtained with different inoculation systems should not be compared. Valid comparisons can only be made among results obtained using the same technique on the same type of material.

## Disposal

- Dispose of contents according to national and local law.
- All patient samples and used components should be handled and disposed of as potentially infectious material.
- Materials of the components:
  - Paper: Instructions for use, patient labels
  - Cardboard: Kit box
  - Plastic: Tubes, caps and dipslides
- When used in accordance with Good Laboratory Practice, good occupational hygiene and the instructions for use, the reagents supplied should not present a hazard to health.

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